

BEAD RUPTOR 24 1.5 ML MICROTUBE CARRIAGE KIT

Sku: 19-010-310S



GENERAL INFO

Bead Ruptor 24 1.5 mL Microtube Carriage Kit

Microtube carriage holds 24 x 1.5 mL microtubes.

Ordering Information: Includes one 24 x 1.5 mL microtube carriage and one finger plate assembly.

SPECIFICATIONS

Capacity: 24 x 1.5 mL microtubes

Volume Range: 200 μ L-1.0mL

Compatibility: most 1.5 mL microtubes