



OMNI
International
The Homogenizer Company™



BEAD RUPTOR 24 2 ML TUBE CARRIAGE KIT

Sku: 19-010-310

GENERAL INFO

Bead Ruptor 24 2 mL Tube Carriage Kit

2 mL tube carriage holds 24 x 2 mL screw cap skirted tubes

Ordering Information: Includes one 24 x 2 mL screw cap tube carriage and one finger plate assembly.

SPECIFICATIONS

Capacity:	24 x 2.0 mL skirted screw cap tube
Volume Range:	200 μ L-1.5mL
Compatibility:	Requires Omni 2.0 mL tubes