

BEAD RUPTOR 24 15 ML TUBE CARRIAGE KIT

Sku: 19-345-015



GENERAL INFO

Bead Ruptor 24 15 mL Tube Carriage Kit

15 mL tube carriage holds 3 x 15 mL conical centrifuge tubes

Ordering Information: Includes one 3 x 15 mL conical centrifuge tube carriage and lock ring.

SPECIFICATIONS

Capacity:	3 x 15 mL conical centrifuge tube
Volume Range:	500 μ L-10mL
Compatibility:	Requires Omni 15 mL tubes
Force Limit:	Maximum force of 6 m/s