

OMNI BEAD RUPTOR 30 ML TUBE CARRIAGE KIT

Sku: 19-345-035



GENERAL INFO

Bead Ruptor 24 30 mL Tube Carriage Kit

30 mL tube carriage holds 6 x 30 mL centrifuge tubes

Ordering Information: Includes one 6 x 30 mL centrifuge tube carriage and one lock ring.

SPECIFICATIONS

Capacity:	6 x 30 mL screw cap conical centrifuge tubes
Volume Range:	5 mL - 30 mL
Compatibility:	Requires Omni 30 mL screw cap tubes
Force Limit:	Maximum force of 5 m/s