

OMNI BEAD RUPTOR 50 ML TUBE CARRIAGE KIT

Sku: 19-345-050



GENERAL INFO

Bead Ruptor 24 50 mL Tube Carriage Kit

50 mL tube carriage holds 3 x 50 mL centrifuge tubes

Ordering Information: Includes one 3 x 50 mL centrifuge tube carriage and one lock ring.

SPECIFICATIONS

Capacity:	3 x 50 mL screw cap conical centrifuge tubes
Volume Range:	5 mL - 50 mL
Compatibility:	Most 50 mL conical screw cap centrifuge tubes
Force Limit:	Maximum force of 5 m/s