


GE-3560 | 1% AGAROSE TAE W/ ETBR LANDSCAPE GEL, 8+1 WELL



Percentage: 1%

Buffer Type: 1X TAE

Gel Type: Agarose

Staining Dye: EtBr

Gel Dimensions: Landscape, 11.7(w) x 5.6(l)cm

Gel Thickness: 6.0mm

Well Format: 8+1

Well Thickness: 1.0mm

Load Volume: 30 μ l

MultiChannel Compatible: Yes, 8 Channel Compatible

Gels per Box: 5 gels/Box

Equipment Compatibility: RunOne or similar size units/systems

Shelf life: 7 months at room temperature