


GE-4710 | 1% AGAROSE TBE W/ ETBR PORTRAIT GEL, 2(12) WELL



Percentage: 1%

Buffer Type: 1X TBE

Gel Type: Agarose

Staining Dye: EtBr

Gel Dimensions: Portrait, 6.5(w) x 9.5(l)cm

Gel Thickness: 6.0mm

Well Format: 2(12)

Well Thickness: 1.0mm

Load Volume: 11µl

MultiChannel Compatible: No

Gels per Box: 10 gels/Box

Equipment Compatibility: RunOne or similar size units/systems

Shelf life: 7 months at room temperature