


## GE-4620 | 1% AGAROSE TBE W/ ETBR MEDIUM GEL, 24+1 WELL



Percentage: 1%

Buffer Type: 1X TAE

Gel Type: Agarose

Staining Dye: EtBr

Gel Dimensions: Medium, 11.7(w) x 8.9(l) cm

Gel Thickness: 6.0mm

Well Format: 24+1

Well Thickness: 1.0mm

Load Volume: 11µl

MultiChannel Compatible: Yes, 12 Channel Compatible

Gels per Box: 5 gels/Box

Equipment Compatibility: RunOne or similar size units/systems

Shelf life: 7 months at room temperature