

Shaker Plate for TubeSpin® Bioreactor 50 – 600 ml



Stable fixation with holding bars



Shaker Plate for TubeSpin® Bioreactor 50

Features of shaker plate:

- Shaker plate fits the relevant shaker model
- Aluminum
- Guarantees largest possible volume per shaker plate
- Secure fixation of the tubes with holding bars

Type	Product-No.	Version	Capacity	Dimensions	Material	Qty/Bag	Qty/Case
			pcs	l x w x h mm		pcs	pcs
	87611	Kühner plate 32 for TubeSpin® Bioreactor 600	8 x 4	800 x 420 x 139	Alu	1	1
	87612	Kühner plate 20 for TubeSpin® Bioreactor 600	5 x 4	500 x 420 x 139	Alu	1	1
	87613	Kühner plate 16 for TubeSpin® Bioreactor 600	4 x 4	420 x 420 x 139	Alu	1	1
	87631	Infors Multitron plate 32 for TubeSpin® Bioreactor 600	8 x 4	850 x 470 x 139	Alu	1	1
	87633	Infors Minitron plate 16 for TubeSpin® Bioreactor 600	4 x 4	480 x 420 x 139	Alu	1	1
	87411	Kühner plate 32 for TubeSpin® Bioreactor 450	8 x 4	800 x 420 x 139	Alu	1	1
	87412	Kühner plate 20 for TubeSpin® Bioreactor 450	5 x 4	500 x 420 x 139	Alu	1	1
	87413	Kühner plate 16 for TubeSpin® Bioreactor 450	4 x 4	420 x 420 x 139	Alu	1	1
	87414	Kühner plate 162 for TubeSpin® Bioreactor 50	18 x 9	800 x 420 x 139	Alu	1	1
	87431	Infors Multitron plate 32 for TubeSpin® Bioreactor 450	8 x 4	850 x 470 x 139	Alu	1	1
	87433	Infors Minitron plate 16 for TubeSpin® Bioreactor 450	4 x 4	480 x 420 x 139	Alu	1	1



1. Shaker plate for 16 x Bioreactor 600
2. Shaker plate for 20 x Bioreactor 600
3. Shaker plate for 32 x Bioreactor 600