

EXPELLPLUS

PIPETTE TIPS

Product Data Sheet



Product Description

| | |
|---------------------|---|
| Catalogue number | 5030042C |
| Product description | ExpellPlus 10/20 µL XL, pre-sterile, clear, bag, individually wrapped |
| Sales Unit | Case / 5 x 400 pcs. |

Product Characteristics

| | |
|-----------------------|---|
| Graduation | Yes |
| Filter | No |
| Filter specifications | N/A |
| Low retention | Yes |
| Sterilised | Yes |
| Autoclavable | Yes, 121 °C, 20 min |
| Purity | Free of RNase, DNase, Human DNA, Pyrogens, PCR Inhibitors |

Size

| | |
|-------------------|----------|
| Max. volume | 20 µL |
| Length of product | 45.80 mm |

Material / Colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | Clear |

Packaging

| | |
|--------------------------|---------------------------|
| Type of packaging | Bag, individually wrapped |
| Pieces/rackbox | |
| Pieces / bag | 400 tips |
| Pieces / pack | |
| Pieces / case | 2000 tips |
| Dimensions case / carton | 232 x 635 x 272 mm |

Storage and Handling

| | |
|-------------|-----------------------|
| Temperature | Below 35°C |
| UV light | Protect from UV light |

Note: *Improper storage can initiate degradation, which results in colour changes and can have negative effects on the physical properties of the products.*

Safety: *The product is not classified as dangerous*

Recycling: *The product is recyclable.*