

Description

- DURAN® GL 45 Multiport Connector Screw Cap, PBT, red

Part number

- 096364

Pieces by sales : 1



Description

- DURAN® GL 45 Multiport Connector Screw Cap, PBT, red
- With 3 tubing connectors, for peristaltic tubing
- Caps have a GL 45 thread that fit DURAN® bottles
- Barbed male hose connector (4.2 mm outer diameter) are suitable for attaching flexible peristaltic PBT tubing with a 3.1 or 3.2 mm 16# inner diameter
- Autoclavable at 121° C or 134 °C
- Construction materials are certified as approved for food contact
- Reusable with a long service life
- Viton seal forms a tight seal with the neck of the bottle
- Thread : GL 45
- Diameter : 54 mm
- Height : 43.25 mm

