


Calibration Statement

Description: VWR Thermometer
Item Number: W0020-110C

This Thermometer has been calibrated and tested against standards traceable to the **National Institute of Standards and Technology (NIST)** using the methods outlined in Monograph 150, Monograph 250-23 from NIST and ASTM E77-07 in manufacture's ISO/IEC and ILAC accredited calibration laboratory.

Each thermometer has a unique serial number for identification.



Andy Franks
Contract Manufacturing Quality Manager
Corning Life Sciences
FranksAJ@Corning.com