

Technical Datasheet

MICROS – MS200 Disposable Microtome Blade

Universal blade for daily routine in histological laboratories, ideally suited for all types of tissue. Short trimming time, very good cutting quality. Especially good for ribbon sectioning cuts, even with very thin cuts. Also suitable for industrial research applications for cutting paraffin and plastic embedded samples. All-rounder with unique hardened cutting edge and excellent durability.



	MS200
Blade Model	MS200 Low profile
Item Number	109211
Measurements (W x D x H)	80 x 8 x 0,25 mm
Blade Angle	35°
Material	stainless steel
Coating	Plasma*
Package size	116 x 53 x 25 mm
Blades / box	50 blades
Blades / sample box	5 blades
CE IVD	✓
Application for	soft and hard tissue, ribbon sectioning, thin sectioning

**By applying plasma and implanting iron along the blade's edge, its hardness can be increased up to 1100 HV, enhancing cutting quality and durability.*

Manual instructions can be found via the QR code.



Revised: 08.01.2025/V1

MICROS Produktions- und Handelsges.m.b.H • Hunnenbrunn - Gewerbezone 5 • A-9300
St. Veit/Glan Telefon +43 4212 30901 • Fax +43 4212 30901-90 • E-Mail office@micros.at

Bank: Unicredit Bank Austria, IBAN: AT211100 0032 7004 6000, BIC: BKAUATWW
Firmenbuch-Nr.: FN75375x, Gerichtsstand: St. Veit/Glan, UID: ATU15049007, EORI: ATEOS1000002137

