

Sterilin Tubes Borosilicate Glass, Culture

Sterilin Culture Tubes are manufactured from Pyrex borosilicate glass, featuring a screw neck



Product Highlights

- Autoclavable

Details

- High resistance to temperature and chemical attack
- Ideal for use with samples sensitive to leaching from plastic tubes
- Tube profile screw thread, culture
- Supplied without caps; autoclavable caps available separately
- Autoclavable at 121°C
- Supplied non-sterile



Sterilin Tubes Borosilicate Glass, Culture

Description	Capacity	O.D. x H	Sterility	Inner Pack Qty.	Case Qty.	Cat. No.
Round bottom tube	7.5mL	13 x 100mm	NS	250	1,000	99449-13
Round bottom tube	11.5mL	16 x 100mm	NS	250	1,000	99449-16
Round bottom tube	15.0mL	16 x 125mm	NS	250	1,000	99449-16X
Round bottom tube	19.0mL	16 x 150mm	NS	250	1,000	99449-16XX
Round bottom tube	24.0mL	20 x 125mm	NS	250	500	99449-20
Round bottom tube	30.0mL	20 x 150mm	NS	250	500	99449-20X
Flat base tube	17.0mL	16 x 125mm	NS	250	1,000	99448-16
Flat base tube	29.5mL	20 x 145mm	NS	250	500	99448-19

Caps to fit Glass Tubes

Description	Liner	Case Qty.	Cat. No.
13mm screw cap with rubber liner (autoclavable)	Rubber	1,000	99999-13
16mm screw cap with rubber liner (autoclavable)	Rubber	1,000	99999-15
20mm screw cap with rubber liner (autoclavable)	Rubber	1,000	99999-18
13mm screw cap with PTFE liner (autoclavable)	PTFE	288	9998-13
16mm screw cap with PTFE liner (autoclavable)	PTFE	288	9998-15
20mm screw cap with PTFE liner (autoclavable)	PTFE	192	9998-18