




Flat tray - 3 L - white

Ref. : G154221



Plus points

Flat trays Largest and most customisable range on the market.

5 colours

5 volumes

Can be used in a wide range of temperatures.

Dimensions suitable for refrigerators.

Can be used as :

- defrosting trays thanks to its wide dividing grid that separates the ice water from the product being defrosted .
 - preparation container for products that drip ingredients tray,
 - containers for preserving comparison samples.
 - Fermenting and decanting containers.
- Nesting and stacking with lid,

Product features

Height	70 mm
Width	230 mm
Length	347 mm
Capacity	3 L
Material	PEHD
Colour	White
Product bottom internal length	280 mm
Product bottom internal width	165 mm
Weight	0,260 kg
GENCODE	3100411542216
Customs code	39241000

Logistic data

Number of items per pallet	1200
Pallet weight	322 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.