




Flat tray - 5 L - white

Ref. : G154321



Plus points

Flat trays Largest and most customisable range on the market.

5 colours

5 volumes

Can be used in a wide range of temperatures.

Dimensions suitable for refrigerators.

Can be used as :

- defrosting trays thanks to its wide dividing grid that separates the ice water from the product being defrosted .
 - preparation container for products that drip ingredients tray,
 - containers for preserving comparison samples.
 - Fermenting and decanting containers.
- Nesting and stacking with lid,

Product features

Height	75 mm
Width	285 mm
Length	437 mm
Capacity	5 L
Material	PEHD
Colour	White
Product bottom internal length	380 mm
Product bottom internal width	220 mm
Weight	0,460 kg
GENCODE	3100411543213
Customs code	39241000

Logistic data

Number of items per pallet	400
Pallet weight	194 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.