




Flat tray - 10 L - white

Ref. : G154521



Plus points

Flat trays Largest and most customisable range on the market.

5 colours

5 volumes

Can be used in a wide range of temperatures.

Dimensions suitable for refrigerators.

Can be used as :

- defrosting trays thanks to its wide dividing grid that separates the ice water from the product being defrosted .
 - preparation container for products that drip ingredients tray,
 - containers for preserving comparison samples.
 - Fermenting and decanting containers.
- Nesting and stacking with lid,

Product features

Height	80 mm
Width	385 mm
Length	535 mm
Capacity	10 L
Material	PEHD
Colour	White
Product bottom internal length	470 mm
Product bottom internal width	315 mm
Weight	0,826 kg
GENCODE	3100411545217
Customs code	39241000

Logistic data

Number of items per pallet	240
Pallet weight	204 kg
Pallet length	800 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.