

## PRODUCT SPECIFICATION

<b>Product Code:</b>	1516/02D
<b>Product Description:</b>	50ml Media Lab Bottles with Cap and Pouring Ring
<b>Brand:</b>	PYREX
<b>Type of Glass:</b>	Borosilicate glass Type 3.3, which conforms to ISO3585 European and US Pharmacopoeia Type 1 glass
<b>Product Specification:</b>	Printed with Trace Code for downloadable batch certificate available on <a href="http://www.scilabware.com">www.scilabware.com</a> White graduations and marking spot Polypropylene cap and drip-free pouring ring  Bottle height 88mm Body diam. 46mm Thread size GL32
<b>Product Standard(s):</b>	The product is manufactured to comply with ISO4796
<b>Pack Quantity:</b>	10 pieces per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of SciLabware Limited by:

Signed:   
\_\_\_\_\_

D Turner  
Quality Technician

Date: 13<sup>th</sup> Feb 2018