

Product Specification

Catalogue No:	1015/22D
Product Description:	1000ml Low Form Beaker
Brand:	Pyrex
Type of Glass:	Type 3.3 Borosilicate Glass, which complies with ISO3585 European and United States pharmacopoeia Glass Type 1
Product Markings:	Product branding, graduated scale, marking spot and trace-code (Colour white) A batch certificate is available from our website www.scilabware.com
Dimensions:	Overall height = 180mm approx. External diameter = 95mm approx.
Standard:	Conforms to ISO3819 and DIN12331
Pack Quantity:	Quantity 10 pieces per pack

SciLabware Limited confirms that the above mentioned product has been manufactured to meet all relevant standards as detailed by our Quality Management System.



R J Turner
(Quality Manager)