



Airtight box - 2.5 L - white

Ref. : G602220



Plus points

Airtight container.
Keeps flavours in and odours out .
Contents are visible and recognizable.
Flexible, watertight lid.
PP container and LDPE lid.

Product features

Height	75 mm
Width	220 mm
Length	220 mm
Capacity	2,5 L
Material	PP
Colour	White
Product bottom internal length	195 mm
Product bottom internal width	195 mm
Product bottom internal height	70 mm
Weight	0.280 kg
GENCODE	3573678654952
Customs code	39241000

Logistic data

Number of items per pallet	360
Pallet weight	133 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.