

# Cole-Parmer®

## Cole-Parmer® Block Heater, 3 Block, Digital, BH-200D-3 & BH-200D-3-HT

- Choice of 130°C or 200°C maximum temperature
- Three block capacity
- Digital setting and display
- Uniform and stable temperatures



# Cole-Parmer® Block Heater, 3 Block, Digital, BH-200D-3 & BH-200D-3-HT



The same precise heating and control of the two block models but with 50% more capacity of a third block. Excellent temperature stability and uniformity is maintained via sensitive electronic controls.

All Cole-Parmer® three block heaters are compatible with the Cole-Parmer® Sample concentrator, for faster sample concentration.

## Technical Specification

Specification	BH-200D-3	BH-200D-3-HT
Number of blocks	3	3
Temperature range	Ambient +8°C to 130°C	Ambient +8°C to 200°C
Temperature stability at 37°C	±0.1°C	±0.1°C
Uniformity within block at 37°C	±0.1°C	±0.1°C
Uniformity within block at 130°C	±1°C	±1°C
Display resolution	0.1°C	0.1°C
Dimensions, mm (w x d x h)	310 x 280 x 115	310 x 280 x 115
Net weight, kg (without blocks)	3.2	3.2
Heater power	450W	450W
Electrical supply	230V, 50Hz	230V, 50Hz
IP Rating	31	31

## Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Block heater, 3 block, digital, 130°C	36610-13	BH-200	BH-200D-3	SBH130D/3
Block heater, 3 block, digital, 200°C	36610-23	BH-200	BH-200D-3-HT	SBH200D/3
Block heater, 3 blocks, digital, 130°C, 120V	36610-12	BH-200	BH-200D-3-120V	SBH130D/3/120V
Block heater, 3 blocks, digital, 200°C, 120V	36610-22	BH-200	BH-200D-3-HT-120V	SBH200D3/120V60

coleparmer.com

USA:  
+1.800.323.4340  
+1.847.549.7600

Canada: +1.800.363.5900  
China: +86.21.5109.9909  
France: +33 (0) 1486 37800  
Germany: +49 (0) 9377 92030

India: +1.800.266.1244  
Italy: +39 (0) 1313 89513  
UK: +44 (0) 1480 272279  
All other countries: +1.847.549.7600

**Cole-Parmer®**  
essentials

an antylia scientific company

Find out more at [antylia.com](http://antylia.com)