



Reinforced container - 10 L - white

Ref. : G110021



Plus points

Reinforced container.
Excellent mechanical resistance: made from low density polyethylene for increased shock resistance.

Product features

Height	160 mm
Width	310 mm
Length	365 mm
Capacity	10 L
Material	PEBD
Colour	White
Product bottom internal length	285 mm
Product bottom internal width	230 mm
Product bottom internal height	150 mm
Weight	0,700 kg
GENCODE	3100411100218
Customs code	39231000

Logistic data

Number of items per pallet	400
Pallet weight	310 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.