



Reinforced container - 19 L - white

Ref. : G109421



Plus points

Reinforced container.
Excellent mechanical resistance: made from low density polyethylene for increased shock resistance.

Product features

Height	170 mm
Width	370 mm
Length	450 mm
Capacity	19 L
Material	PEBD
Colour	White
Product bottom internal length	365 mm
Product bottom internal width	285 mm
Product bottom internal height	165 mm
Weight	1,700 kg
GENCODE	3100411094210
Customs code	39231000

Logistic data

Number of items per pallet	240
Pallet weight	278 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.