



Reinforced container - 30 L - white

Ref. : G113721



Plus points

Reinforced container.
Excellent mechanical resistance: made from low density polyethylene for increased shock resistance.

Product features

Height	200 mm
Width	440 mm
Length	540 mm
Capacity	30 L
Material	PEBD
Colour	White
Product bottom internal length	450 mm
Product bottom internal width	340 mm
Product bottom internal height	190 mm
Weight	1,670 kg
GENCODE	3100411137214
Customs code	39231000

Logistic data

Number of items per pallet	160
Pallet weight	298 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.