



## Reinforced container - 45 L - white

Ref. : G155021



### Plus points

Reinforced container.  
Excellent mechanical resistance: made from low density polyethylene for increased shock resistance.

### Product features

Height	305 mm
Width	445 mm
Length	550 mm
Capacity	45 L
Material	PEBD
Colour	White
Product bottom internal length	415 mm
Product bottom internal width	305 mm
Product bottom internal height	295 mm
Weight	2.640 kg
GENCODE	3100411550211
Customs code	39231000

### Logistic data

Number of items per pallet	120
Pallet weight	316 kg
Pallet length	800 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.