



Reinforced container - 55 L - white

Ref. : G155821



Plus points

Reinforced container.
Excellent mechanical resistance: made from low density polyethylene for increased shock resistance.

Product features

Height	380 mm
Width	430 mm
Length	555 mm
Capacity	55 L
Material	PEHD
Colour	White
Product bottom internal length	445 mm
Product bottom internal width	320 mm
Product bottom internal height	370 mm
Weight	3,000 kg
GENCODE	3100411558217
Customs code	39231000

Logistic data

Number of items per pallet	48
Pallet weight	147 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.