



Square basin - 1 L - white

Ref. : G605621



Plus points

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect .
Moulded rim makes it easy to grasp.

Product features

Height	65 mm
Width	160 mm
Length	160 mm
Capacity	1 L
Material	PEBD
Colour	White
Product bottom internal length	120 mm
Product bottom internal width	1202 mm
Weight	0,090 kg
GENCODE	3100416056213
Customs code	39241000

Logistic data

Number of items per pallet	2160
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.