



Square basin - 3 L - white

Ref. : G605821



Plus points

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect .
Moulded rim makes it easy to grasp.

Product features

Height	90 mm
Width	220 mm
Length	220 mm
Capacity	3 L
Material	PEBD
Colour	White
Product bottom internal length	160 mm
Product bottom internal width	160 mm
Weight	0,210 kg
GENCODE	3100416058217
Customs code	39241000

Logistic data

Number of items per pallet	1200
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.