



Square basin - 4.5 L - white

Ref. : G605921



Plus points

- Basin fit for food contact.
- Highly durable and chemical resistant.
- Easy to clean.
- Nesting with no suction effect .
- Moulded rim makes it easy to grasp.

Product features

Height	110 mm
Width	260 mm
Length	260 mm
Capacity	4,5 L
Material	PEBD
Colour	White
Product bottom internal length	195 mm
Product bottom internal width	200 mm
Weight	0,280 kg
GENCODE	3100416059214
Customs code	39241000

Logistic data

Number of items per pallet	640
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.