



Square basin - 7 L - white

Ref. : G606021



Plus points

Basin fit for food contact.
Highly durable and chemical resistant.
Easy to clean.
Nesting with no suction effect .
Moulded rim makes it easy to grasp.

Product features

Height	125 mm
Width	300 mm
Length	300 mm
Capacity	7 L
Material	PEBD
Colour	White
Product bottom internal length	225 mm
Product bottom internal width	225 mm
Weight	0,490 kg
GENCODE	3100416060210
Customs code	39241000

Logistic data

Number of items per pallet	540
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.