



Square basin - 10 L - white

Ref. : G606121



Plus points

- Basin fit for food contact.
- Highly durable and chemical resistant.
- Easy to clean.
- Nesting with no suction effect .
- Moulded rim makes it easy to grasp.

Product features

Height	130 mm
Width	350 mm
Length	350 mm
Capacity	10 L
Material	PEBD
Colour	White
Product bottom internal length	270 mm
Product bottom internal width	270 mm
Weight	0,620 kg
GENCODE	3100416061217
Customs code	39241000

Logistic data

Number of items per pallet	360
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.