



## Square basin - 12 L - white

Ref. : G606221



### Plus points

Basin fit for food contact.  
Highly durable and chemical resistant.  
Easy to clean.  
Nesting with no suction effect .  
Moulded rim makes it easy to grasp.

### Product features

Height	140 mm
Width	390 mm
Length	390 mm
Capacity	12 L
Material	PEBD
Colour	White
Product bottom internal length	290 mm
Product bottom internal width	290 mm
Weight	0,710 kg
GENCODE	3100416062214
Customs code	39241000

### Logistic data

Number of items per pallet	240
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.