



## Round basin - 7.5 L - white

Ref. : G606721



### Plus points

- Basin fit for food contact.
- Highly durable and chemical resistant.
- Easy to clean.
- Nesting with no suction effect .
- Moulded rim makes it easy to grasp.

### Product features

Height	125 mm
Diameter	350 mm
Capacity	7,5 L
Material	PEBD
Colour	White
Product internal diameter	255 mm
Weight	0.460 kg
GENCODE	3100416067219
Customs code	39241000

### Logistic data

Number of items per pallet	540
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.