



Round basin - 14 L - white

Ref. : G606921



Plus points

Round dividing dough container.
Highly resistant to changes in temperature.
Highly flexible and shock-resistant.
Compatible with tray holders.

Product features

Height	175 mm
Diameter	430 mm
Capacity	14 L
Material	PEBD
Colour	White
Product bottom internal height	170 mm
Product internal diameter	360 mm
Weight	0,730 kg
GENCODE	3100416069213
Customs code	39241000

Logistic data

Number of pieces by subpacking unit	20
Subpacking weight	14,6 kg
Subpacking length	445 mm
Subpacking width	415 mm
Subpacking height	400 mm
Number of items per pallet	360
Pallet weight	272 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	1390 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.