



## Pedal container with yellow lid - 100 L

Ref. : G614425



### Plus points

- Pedal container.
- Equipped with two wheels.
- Easy to move around.
- Hands-free operation for sanitary disposal.
- Easy to clean.
- Coloured lids for selective waste sorting.

### Product features

Height	922 mm
Width	415 mm
Length	472 mm
Capacity	100 L
Material	PP
Colour	Yellow
Weight	8.000 kg
GENCODE	3573678648579
Customs code	39241000

### Logistic data

Number of items per pallet	8
Pallet weight	74 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	2050 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.