



Deep rectangular 16 L basin - white

Ref. : G606321



Plus points

- Basin fit for food contact.
- Highly durable and chemical resistant.
- Easy to clean.
- Nesting with no suction effect .
- Moulded rim makes it easy to grasp.

Product features

Height	200 mm
Width	350 mm
Length	400 mm
Capacity	16 L
Material	PEBD
Colour	White
Product bottom internal length	290 mm
Product bottom internal width	280 mm
Weight	0,950 kg
GENCODE	3100416063211
Customs code	39241000

Logistic data

Number of items per pallet	270
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.