

Petri dish

Reference	Ø x height (mm)	Pack unit
5058501	12 x 40	40
5058541	15 x 60	30
5058543	15 x 80	36
5058546	15 x 90	36
5058544	10 x 100	36
5058545	15 x 100	36
5058548	20 x 100	36
5058550	20 x 120	72
5058558	25 x 150	36
5058559	30 x 180	6
5058560	30 x 200	6
5058561	50 x 200	6



SODA LIME

Petri dish borosilicate 3.3

Reference	Ø x height (mm)	Pack unit
5058641	15 x 60	28
5058643	15 x 80	36
5058646	15 x 90	36
5058645	15 x 100	36
5058648	20 x 100	30
5058650	20 x 120	68
5058658	25 x 150	36
5058659	30 x 180	6
5058660	30 x 200	5
5058661	50 x 200	5



121°C

BORO 3.3

Pasteur pipette

Reference	Outside Ø (mm)	Length (mm)	Tip	Pack unit
5426015	6,5-8	145-150	open	250
5426023	6,5-8	225-230	open	250
5426027	6,5-8	270	open	250
5426015F	6,5-8	145-150	closed	250
5426023F	6,5-8	225-230	closed	250
5426027F	6,5-8	270	closed	250



SODA LIME

Pasteur pipette cotton plugged

Reference	Outside Ø (mm)	Length (mm)	Tip	Pack unit
5426015C	6,5-8	145-150	open	250
5426023C	6,5-8	225-230	open	250
5426027C	6,5-8	270	open	250
5426015CC	6,5-8	145-150	closed	250
5426023CC	6,5-8	225-230	closed	250
5426027CC	6,5-8	270	closed	250



SODA LIME