

## TECHNICAL DATASHEET

### Petri dish

Petri dish is a shallow, cylindrical, round glass dish, which is used to culture different microorganism and cells  
Base and lid are flat inside and out, and free from bubbles and streaks for viewing under a microscope to be sharp

Cat nr	Ø x height (mm)
5058501	12x40
5058541	15x60
5058543	15x80
5058546	15x90
5058544	10x100
5058545	15x100
5058548	20x100
5058550	20x120
5058558	25x150
5058559	30x180
5058560	30x200
5058561	50x200



#### Material:

Clear soda-lime glass

#### Features:

Microorganism and cell culture  
Hydrolytic resistance class III  
Medium thermal shock endurance