

# Cole-Parmer®

## Cole-Parmer® Colony Counter CC-200

- Pressure sensitive counting
- Average count facility
- Bright white energy saving LED lighting
- With BioCote antimicrobial protection
- Audible count confirmation
- Choice of light or dark background
- Connectivity to printer or computer



# Cole-Parmer® Colony Counter CC-200



Touch pressure with felt tip marker on petri dish registers cumulative count on the digital display with confirmation by audible tone (can be turned on or off). The pressure required to register a count can be adjusted to suit each user.

Averaging facility calculates average count over multiple plates. Counting results as well as useful statistics including SD can be sent directly to the accessory printer or to a computer via a USB cable supplied.

Sub-stage illumination by low energy bright LED's allows glare-free optimum viewing. A switchable black background is provided to enhance viewing of translucent and difficult to see colonies. Supplied with two Wolffhuegel graticules and dish centering adaptors to facilitate use with 50mm to 90mm dishes. A choice of magnifiers and a printer are available as optional accessories.

## Technical Specification

| Specification          | CC-200                    |
|------------------------|---------------------------|
| Lighting               | White LED array           |
| Digital display        | 3 digit LED               |
| Count                  | 0 to 999                  |
| Dimensions (w x d x h) | 310 x 300 x 140mm         |
| Net weight             | 1.5kg                     |
| Electrical supply      | 90 to 240V, 50/ 60Hz, 70W |

## Ordering Information

| Description   | Ordering Number | Series No. | Model No. | Legacy Sku.  |
|---|-----------------|------------|-----------|--------------|
| Colony counter, advanced                                | 14212-04        | CC-200     | CC-200    | SC6PLUS      |
| Colony counter, advanced with IQ/OQ documentation       | 14212-06        | CC-200     | CC-200    | SC6PLUS/IQOQ |
| 1.7x magnifier  | 14212-60        | CC-200     | N/A       | SC6/ 1       |
| 3x magnifier  | 14212-61        | CC-200     | N/A       | SC6/ 1/ 3    |
| Wolffhuegel graticule & segmentation discs (pack of 10) | 99952-84        | CC-200     | N/A       | SC6/ 2       |
| Spare dish centering adaptors (pack of 2)               | 99952-85        | CC-200     | N/A       | SC6/ 3       |
| Clear protective discs (pack of 5)                      | 99952-86        | CC-200     | N/A       | SC6/ 4       |
| Spare receiving dish                                    | 99952-87        | CC-200     | N/A       | SC6/ 5       |
| Accessory Printer                                       | 03012-99        | CC-200     | N/A       | SMP30/ 1     |

coleparmer.com

USA:  
+1.800.323.4340  
+1.847.549.7600

Canada: +1.800.363.5900  
China: +86.21.5109.9909  
France: +33 (0) 1486 37800  
Germany: +49 (0) 9377 92030

India: +1.800.266.1244  
Italy: +39 (0) 1313 89513  
UK: +44 (0) 1480 272279  
All other countries: +1.847.549.7600

**Cole-Parmer®**  
essentials

an antylia scientific company

Find out more at antylia.com