


## Product Specification

<b>Product Code:</b>	2115/02M
<b>Product Description:</b>	Glass Weighing Scoop
<b>Brand:</b>	PYREX
<b>Glass Type:</b>	Borosilicate Glass Type 3.3 which complies with ISO3585 European and United States Pharmacopoeia Glass Type 1
<b>Product Specification:</b>	Capacity – 3ml Overall length – 70mm
<b>Pack Quantity:</b>	10 pieces per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001:2008.

Signed on behalf of SciLabware Limited by:



Signed: \_\_\_\_\_  
K. Skupniewicz  
Quality Engineer

Date: 23 Aug 2017