

## Product Specification

Product Code	CP0010P
Product Description	AZLON MEASURING CYLINDER, CLASS B, MOULDED GRADUATIONS, PP, 10 ML
Material particular	Polypropylene
Capacity ml	10
Colour	Ultra Clear PP
Autoclavable	Yes, See <a href="http://www.scilabware.com">www.scilabware.com</a>
Pack quantity	10

## Dimensions and tolerances

Height mm	145
Diameter mm	14
Graduation ml	20.0

## Standards

Conforms to	ISO 6706
-------------	----------

## Product Compliance

TSE/BSE Compliant	Yes
Reach	Yes
Latex	Yes

**SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO 9001:2008**

**Edel Tuite**

**Azlon Quality Manager**

Date issued 1/19/2017