


Product Specification

Product Code:	5810/3
Product Description:	Base for Filtration System, 47mm Diameter Membranes
Brand:	PYREX
Glass Type:	Borosilicate Glass Type 3.3 which complies with ISO3585 European and United States Pharmacopoeia Glass Type 1
Product Markings:	Product branding, country of origin, product code and socket size.
Product Specification:	Ground Joint (Socket Size) - 40/38 Sintered Disc Porosity Grade 3 to support the filter membrane Fully autoclavable
Pack Quantity:	1 piece per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001:2008.

Signed on behalf of SciLabware Limited by:



Signed: _____
K. Skupniewicz
Quality Engineer

Date: 23 Aug 2017