


## Product Specification

<b>Product Code:</b>	5810/5
<b>Product Description:</b>	Clip for Filtration System, 47mm Diameter Membranes
<b>Brand:</b>	PYREX
<b>Product Specification:</b>	Metal retaining clip Fully autoclavable
<b>Pack Quantity:</b>	1 piece per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001:2008.

Signed on behalf of SciLabware Limited by:



Signed: \_\_\_\_\_  
K. Skupniewicz  
Quality Engineer

Date: 23 Aug 2017