

## EU Declaration of Conformity

<b>Product</b>	Laboratory Equipment	<b>File Number</b>	TFR909081_001
<b>Manufacturer</b>	Cole-Parmer Ltd Beacon Road Stone, Staffordshire ST15 0SA United Kingdom		TFR909081_002

**This declaration of conformity is issued under the sole responsibility of the manufacturer**

**Object of Declaration**                      Shaker Incubator  
*(reference the attached list of catalogue numbers)*

*The object of the declaration described above is in conformity with the relevant Union Harmonisation Legislation:*

<b>Low Voltage Directive</b>	2014/35/EU
<b>EMC Directive</b>	2004/108/EC
<b>RoHS Directive</b>	2011/65/EC

*References to the relevant harmonised standards used or references to the other technical specifications in relation to which conformity is declared:*

IEC/EN 61010-1:2010	Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1: General requirements.
IEC/EN 61010-2-010:2003	Particular requirements for laboratory equipment for the heating of materials.
IEC/EN 61010-2-051:2015	Particular requirements for laboratory equipment for mixing and stirring.
IEC/EN 61326-1:2006	Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 1: General requirements (Class A).

**Signed for and on behalf of the above manufacturer**

<b>Additional Information</b>	Year of CE Marking: 2013
<b>Place of Issue</b>	Stone, Staffordshire, UK
<b>Date of Issue</b>	29 April 2013. Revised 09 March 2017
<b>Authorised Representative</b>	Steve Marriott
<b>Title</b>	Technical Director
<b>Signature</b>	

Catalogue Number	Model	Description
79520-00	SI-200D-120	Shaking Incubator, 51 L capacity, 120 V, 60 Hz
79520-04	SI-200D	Shaking Incubator, 51 L capacity, 230 V, 50 Hz
79520-08	SIM-200D-M-120	Microplate Shaking Incubator, 120 V, 60 Hz
79520-10	SIM-200D-M	Microplate Shaking Incubator, 230 V, 50 Hz
79520-14	SI-200D-L-120	Large Capacity Shaking Incubator, 115 L capacity, 120 V, 60 Hz
79520-16	SI-200D-L	Large Capacity Shaking Incubator, 115 L capacity, 230 V, 50 Hz
79520-20	SIC-200D-C-120	Large Capacity Refrigerated Shaking Incubator, 115 L capacity, 120 V, 60 Hz
79520-22	SIC-200D-C	Large Capacity Refrigerated Shaking Incubator, 115 L capacity, 230 V, 50 Hz