

Corning Incorporated
Life Sciences

Registered
ISO 9001:2008

Product Description

Catalog Number: 2593

Product Description: Costar® 1 x 8 Stripwell™, medium binding EIA/RIA plate, flat bottom, without lid

Component Materials:

- Strips - Virgin Polystyrene, meets *USP, Class VI* requirements for plastic containers and closures.
- Holder - Polystyrene, Heavy metal free (meets *CONEG* req.) color concentrate.

Product Dimensions:

Length of Plate	-	5.025 in.	Diameter of Well @ top	-	.270 in.
Width of Plate	-	3.370 in.	Diameter of Well @ bottom	-	.250 in.
Depth of Well	-	.420 in.	Height	-	.560 in.
Tolerances of Dimensions	-	+/- .010 in.	Working Volume per well	-	.36 mL

Surface Characterization:

Surface is characterized to be medium binding for EIA/RIA tests. Binding capacity is approximately 250 ng IgG/cm².

Optical Characteristics:

The product is made of clear polystyrene to allow background absorbance to be $\pm .005$ optical density units from the mean.

Binding Characteristics:

Each lot has been tested for the attribute of well to well protein binding utilizing a modified ELISA to measure precision and percent low and high wells from the mean. Coefficient of variation must be $\leq 5\%$ and high and low wells $\leq 15\%$ from the mean. Product has a Costar Certified Surface Chemistry.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

- Visual Attributes: Visual and microscopic examination of the product.
- Packaging: Inspection for accurate labeling, and correct product configuration.
- Functional: Tested for protein binding precision.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 4