

Kimtech™ Science Wipes



Low-lint, **high absorbency** and anti-static

Chemically inert wipes for delicate tasks

Convenient pop-up box to ensure **clean and fast dispensing**

Kimtech™ Science Wipes are designed to maintain cleanliness in your laboratory or working environment. Effective cleaning of equipment, materials and working environments is important to ensure the accuracy, stability and longevity of your lab. Kimtech™ Science provides wipes for delicate tasks.

The light duty wipes are suitable for a variety of wiping tasks such as the cleaning of delicate glassware, equipment and surfaces. Lenses, slides and optical instruments can be wiped with ease, and the wipes are suitable for ISO Class 4 or higher laboratory environments.

Made from 100% pure fibre, giving high absorbency and a soft feel, the white-coloured wipes are chemically inert allowing use with solvents and chemicals. The wipes are also dispensed from a flexible pop-up box, to ensure clean and fast dispensing, even while on the move. Suitable for use within laboratories in manufacturing and healthcare, as well as industry sectors including optical, laser and electronics, the wipes feature low-lint tissue softness, making them ideal for delicate cleaning applications. They are also suitable for cleaning liquids and dust, and are equipped with anti-static dispensing to reduce lint and antistatic discharge.



Key Features

- For delicate equipment, surfaces, optics and materials, our low-lint and high absorbency white wipes have you covered
- Light duty, white-coloured wipes made from 100% pure fibre
- Chemically inert and anti-static
- Provided in a hygienic, user-friendly pop-up box
- Suitable for use in ISO Class 4 or higher environments
- Suitable for applications in healthcare, manufacturing, optics and electronics

Product Specifications

Kimtech™ Science Delicate Task Wipes				
Code	Ply	Size (WxL)	Colour	Total/Case
7557	2-ply	20 cm x 20.5 cm	White	24 boxes x 100 sheets = 2,400 sheets
7558	2-ply	20 cm x 20.5 cm	White	15 boxes x 200 sheets = 3,000 sheets

Kimtech™ Science Precision Wipes				
Code	Ply	Size (WxL)	Colour	Total/Case
7551	1-ply	28.5 cm x 31 cm	White	15 boxes x 198 sheets = 2,970 sheets
7552	1-ply	11.2 cm x 20.8 cm	White	30 boxes x 286 sheets = 8,580 sheets