

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

## Product Description

Catalog Number: 3008

Product Description: Costar® Cell Lifter

**Component Materials:**

Lifter - High density polyethylene, meets *USP, Class VI* requirements for plastic containers and closures.

**Product Dimensions:**

Length of handle	- 7.36 in.	Length of blade	- 0.75 in.
Tolerances of Dimensions	- +/- 0.05 in.		

**Sterilization:**

The lot has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137 (TIR 33) *Sterilization of healthcare products-Requirements for validation and routine control-Radiation sterilization*.

Sterility Assurance Level: SAL 10<sup>-3</sup>

**Pyrogens:**

The product has been tested and has met the criteria established in ANSI/AAMI ST 72:2002/(R) 2010: *Bacterial endotoxins - Test methodologies, routine testing, and alternative to batch testing*. Results: less than 0.1 EU/mL

**Performance Testing:**

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes: Visual examination of the product.

Packaging: Inspection for seal and barrier integrity, accurate labeling, and correct product configuration.

**Lot Number Designation:**

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification. The entire lot is sterilized as one unit.

Rev No: 6