


## FICHE TECHNIQUE

---

<b>Référence</b>	003769
<b>Désignation</b>	Corning® 384-well Flat Clear Bottom



---

### Description :

- ✓ Plate Format 384-well
- ✓ Plate Feature Low flange
- ✓ Plate Color Clear bottom white
- ✓ Well Bottom Flat
- ✓ Well Bottom Color Clear
- ✓ Well Volume 95  $\mu\text{L}$
- ✓ Cell Growth Area 0,06  $\text{cm}^2$
- ✓ Recommended Medium Well Volume 25 - 50  $\mu\text{L}$
- ✓ Recommended Working Volume 20-80  $\mu\text{L}$
- ✓ Sterile Yes
- ✓ Lids Included Oui
- ✓ Surface Coating Poly-D-Lysine

### Unité de vente :

100 unités

---