

**Corning Incorporated
Life Sciences**

Registered
ISO 9001:2008

Product Description

Catalog Number: 3895

Product Description: Corning® 1536-well, non-binding, black plate with clear flat bottom, without lid

Component Materials:

- Plate - Virgin Polystyrene, meets *USP, Class VI* requirements for plastic containers and closures. Black concentrate
- Plate bottom - Virgin Polystyrene, meets *USP, Class VI* requirements for plastic containers and closures.

Product Dimensions:

- | | | | | | |
|--|---|--------------|-------------------------------|---|----------------------------------|
| Length of Plate | - | 5.030 in. | Width of Square Well @ Top | - | .067 in. |
| Width of Plate | - | 3.365 in. | Width of Square Well @ Bottom | - | .0035 in. / .022 cm ² |
| Depth of Well | - | .195 in. | Plate Height | - | .410 in. |
| Tolerances of Dimensions | - | +/- .010 in. | Total well volume | - | 12.2 µL |
| Recommended working volume per well is | - | up to 10 µL | | | |

Surface Characterization:

Surface is characterized to be non-binding.

Binding Characteristics:

Each lot has been tested for low well to well protein binding utilizing a standard protein binding technique.

Optical Characteristics:

The product is made of opaque black polystyrene walls to minimize well to well crosstalk and background fluorescence and /or luminescence. The bottom is made of clear polystyrene to allow visual clarity for optical density or bottom read fluorescent / luminescent measurements.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Visual Attributes:

Visual examination of the product.

Packaging:

Inspection for accurate labeling, and correct product configuration.

Functional:

Standard protein binding technique.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian Date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 1