


## FICHE TECHNIQUE

---

<b>Référence</b>	004583
<b>Designation</b>	Plaque 384 puits à fond en verre Corning BioCoat™ - traitement collagène I



---

### Description :

Qty./Pk 1 / Pk  
Qty./Cs 10 / Cs  
Brand Corning®  
Plate Format 384-well  
Plate Color Clear bottom black  
Well Bottom Flat  
Well Bottom Color Clear  
Well Volume 82 µL  
Sterile No  
Lids Included Oui  
Storage 2°C- 8°C  
Surface Coating Collagen I